W3 Scan

Sistema Digital de Diagnóstico por Ultrasonido Wireless
Digital Ultrasound Diagnostic Wireless System

IMPORTANTE Mantenimiento Sonda W3 SCAN

1. Recomendaciones de uso

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IMPORTANT W3 Scan probe maintenance

5. Correct use recommendations

- To avoid probe overheating issues, **place the probe on stand-by for at least 5 minutes after prolonged use (25min)**
- Turn the probe off when not in use, **never leave the probe on**
- After performing the scan, **clean the probe carefully**. Use a wet wipe to avoid **scratching the porous membrane**.
- Use original charger when charging the probe.
- Handle the probe with care.

6. Storage conditions

- Temperature: -40 to +50°C
- Humidity: 10 to 85% (without condensation)
- Atmospheric pressure: between 500 and 1060hPa

7. Working conditions

- Temperature: +5 to +40°C
- Humidity: 30 to 85% (without condensation)
- Atmospheric pressure: between 700 and 1060hPa

8. Transport conditions

- Temperature: +5 to 40°C
- Humidity: 30 to 85% (without condensation)
- Atmospheric pressure: between 700 and 1060hPa

1. INTRODUCTION

Conventional probes need to connect to a screen via cables. The W3 SCAN is equipped with the latest in data transfer technology and imaging, not requiring cables for it to function.

The W3 SCAN probe connects without cables via Wi-Fi and an android application (APP). This allows for tables, mobile phones, or any Android device to connect without the need for cables or internet connection. The *Go Scan* APP has been developed for intuitive use with the probe.

We recommend before using the KUBSU W3 probe, to read this manual in detail.

1.3 REPRODUCTIVE ULTRASOUND SCANNING

W3 Scan farm applications

- Early pregnancy diagnostics
- Pubertal vs non-pubertal determination
- Back fat measurement
- Pathology diagnostics

2019/05/31 11:14

Depth: 160mm

Gain: 110 dB

BLADDER

UTERUS

EMBRYONIC
VESICLE


OVARY WITH
LUTEAL BODIES

VENOSUS PLEXUS

INTESTINAL LOOP

Freeze

22 days

1.4 DEVICE CHARACTERISTICS

Image Display Device	Monitor or Smart phone (Android)
Probe Type	80° Sectorial Probe
Working Frequency	3.5 MHz
Scanning depth	100 – 180 mm
Viewing mode	B
Grey Scale	256
Battery Life	> 3 hours
Size	140 x 45 x 40mm
Weight	< 0,2 Kg

WARNING: After 20-25min of continuous use, the probe may overheat, leading to reduction in image quality or break down. We recommend putting the probe on Stand-by between each diagnostic (Scan). A 5min cooldown period every 25min is recommended.

2. BEFORE SCANNING

Please read this manual carefully to avoid unnecessary risks during the use of the device.

2.6 PACKAGE OPENING

To avoid damage during transport, the W3 probe is carefully shipped in a padded case. Before opening the box, please verify that the package has arrive in good conditions, and then revise the contents to ensure nothing is missing.

W3 Contents List

Item	Included
W3 Probe	☐
USB Charger	☐
Wrist strap	☐
Probe Bag	☐
Charger	☐

Even though the W3 probe has been carefully packaged, verify that no item is damage or defects can be observed. If damage or defects are found, contact your distributor.

2.7 APP INSTALATION

If the W3 APP is not installed, contact the local distributor or the KUBUS technical service (+34 616 016 799).

Installation instructions.



5. Download the APP from the link given by the KUBUS technical service.
6. Transfer the APP to Tablet or Smart Phone
7. Open the APP from the Tablet or Smart Phone. Search for the APP named “**Go-Scan**”.



2.8 CONNECT THE APP WITH W3 PROBE

Maintained the probe button pressed, it will automatically connect with the APP. If the connection were to fail at any point, the probe connection should be reconfigured following these simple steps:

1. Open **Configuration** menu of the device (Tablet or Smart Phone)
2. Open the WLAN configuration menu (as shown in the images below)
3. Click on the Probe ID. The ID varies but as an example, the ID will be similar to “UScanYAID031”.

8. The password is the same as the probe's ID (be aware of capital letters). After inputting the password, allow the system to connect.

9. Finally, open the "Scan Go" APP, the probe will be paired with the application.

2.9 CONNECT TO THE APP WITH THE W3 PROBE

Hold down the start button to start-up the probe. After turning the probe on, the Wi-Fi and battery indicators will light up. The battery indicator will show the battery level. After a few seconds, the Wi-Fi probe will turn on, indicating that the probe is ready to work with.


To turn off the probe, hold down the power button. The battery indicator will turn off when the probe switches off.

2.10 WIRELESS CONNECTION

The Wireless Wi-Fi indicator will flash until it finishes connecting. Once connected, the light will be solid, and the probe is ready to use.

3. WORKING WITH THE APP

3.1 Ultrasound scanning

13. KUBUS logo

14. Gain scale bar

15. Freeze button

16. Configuration

17. View data

18. Save image

19. Save video

20. Play/Pause button

21. Image advance bar

22. Colour modification

23. Area of image

24. Date, time, Depth and gain indicator

Until the probe turns on, the screen image will be black. To obtain an image, press the freeze button on the screen, or the central button of the probe. The probe will generate ultrasound images.

The image gain can be modified in the adjustment bar (2, as shown in the image) Image depth can be adjusted, by pinching the image up or down, until the desired image depth is obtained.

When the image is frozen, the display will indicate the number of images taken during the scan. The sequence of images can be viewed pressing the play button (8). To access a specific image within a sequence, whilst the scanner is paused, use the advance bar (9) to find the desired image.

The four left side buttons are configuration (4), view data (5) save image (6) and save video (7)


3.2 Save Image or video

Whilst the screen is paused, press the save image button (6) to generate an image file. This can be accessed by pressing the view data button (5). For video, the same procedure can be done, except that to save the video, the save video button (7) must be pressed.


4. TRANSFER IMAGES AND VIDEOS



Images stored in the device can be sent to another device. This can be done via USB, e-mail or or an app that allows sending images and documents. Image storage can be seen via the device (1) or through the generated files (2). The program can generate Excel and PDF files, as well as send and share them.

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

5. MAINTENANCE

5.1 Charging the probe

Charge the probe when battery is low via its USB charger. Turn off the probe previously by pressing a few seconds the centre button. You can use any USB port for charging (PC, wall, or power bank). During charge, the indicator light will turn red. When the light turns blue, the battery is fully charged. Once the charging process is completed, disconnect the USB, and seal the charge port with the watertight seal to protect from water.

5.2 Cleaning

After each scan, clean the outer layer of the probe with tap water, a wet towel or wet wipes. The probe can be cleaned with running water, but never submerge.


5.3. Storage

When we have finished using the probe, store with care in its bag and in its padded case to avoid bumps and knocks. Store at a temperature between 0-40°C

6. BACK FAT MEASUREMENT

To carry out correct back fat measurements, the gain should be lowered as well as the scanning depth slightly. Freeze the image and with the measuring tool highlight the region of backfat to measure. The thickness of the backfat will show automatically.

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KUBUS LAB S.A.