

KUBUS

QS 1000 PACKAGING MACHINE AND AUTOMATIC BLISTER SEALING MACHINE

REF: 3780

QS1000 Features

- 1000blister/h production capacity.
- Complete process: Cut, fill, seal and label semen doses with a simple touch of the screen.
- Optical-laser sensor to cut and fill semen doses for SAI or PCAI.
- Integrated PLC touchscreen to select fill volume and number of units.
- The body is made from smooth finished fibreglass, providing a resistant and low weight material.
- The body is easy to clean reducing contamination problems with a wet cloth.
- Customisable labels can be placed directly on the blister in the process, where the user can personalise logos and text which appear on the labels.



TECHNICAL SPECIFICATIONS

- Windows 10 Pro integrated PC.
- PLC touchscreen.
- Multilingual.
- Emergency shutdown button.
- Silicone tubes can be exchanged screwlessly in two simple steps.
- USB connectivity or Ethernet remote connectivity.
- Operated with air compressor.
- Power: 250 W
- Voltage: 220V
- Weight: 80 kg
- Size: 800x880x750cm

CONTACT US

 www.kubus-sa.com

 kubus@kubus-sa.com

 +34 91 636 02 68